

Canon

CX-1

Digital Retinal Camera **MYD / NM**

Mydriatic/Non-mydriatic Hybrid Retinal Camera



# A simple, compact unit for high-quality mydriatic and non-mydriatic photography.

The CX-1 provides one-touch operation to select all modes, and features mydriatic and non-mydriatic photography.

Its sleek, compact body boosts efficiency in exam rooms.

It features a digital camera unit combining Canon's SLR camera technology and fundus camera technology. All this delivers high-definition, high-quality retinal images.





## One Unit for Simple Mydriatic/ Non-mydriatic Photography

CX-1 is a single system capable of mydriatic and non-mydriatic photography. Mydriatic photography uses an optical viewfinder, while non-mydriatic photography uses the monitor for fundus observations. Modes can be switched with one touch of the panel button.

## Digital Camera Unit for Fundus Camera

A digital camera unit has been developed specifically for the CX-1, combining the functions of SLR camera technology and fundus camera technology. The result is clear with greater-detail, higher-quality images.

## Mydriatic/Non-mydriatic Color Low-light Photography

Color photography with the CX-1 is possible with low flash intensity for less strain on the patient, regardless of mydriatic and non-mydriatic modes.

## Five Photography Modes

Five photograph modes: COLOR/RED FREE/COBALT/FLUO/FAF. Both mydriatic and non-mydriatic images can be captured. All modes can be selected by simply using the button on the panel.

## Non-mydriatic Fundus Autofluorescence (FAF) Photography

The CX-1 can take fundus autofluorescence (FAF) photography for non-mydriatic observations. Imaging is quick and easy, and also less stressful for patients.





## Intuitive Operation

Lightweight, compact design means the CX-1 is a breeze to handle, and easy to assist the patient in opening their eyelids. With pan and tilt functions, it can be stopped in the direction you want to photograph.


# A simple, compact retinal camera with functions

## Five Photography Modes

COLOR / RED FREE / COBALT / FLUO / FAF

Each mode is available with mydriatic and non-mydriatic photography.

FLUO photography can be used in 1 photo/sec burst mode.



## Easy-to-use Control Panel

Frequently used functions, like switching between mydriatic and non-mydriatic observations or selecting modes, are arranged intuitively. The result is simple and easy-to-use operations.



# for mydriatic and non-mydriatic observations




## Dedicated Digital Camera for Retinal Photography

A digital camera unit has been developed specifically for the CX-1, combining Canon's expertise with SLR camera technology and fundus camera technology. The CX-1 and Canon's digital camera technology operate in concert to deliver high-quality images.



# Five photography modes available in both mydriatic and non-mydriatic modes.

## CX-1 Image samples



## Fundus Autofluorescence (FAF) Photography Available for Both Mydriatic and Non-mydriatic


The CX-1 comes standard with the fundus autofluorescence (FAF) photography function for mydriatic and non-mydriatic observation of both eyes. Results are available quickly and easily, which can also reduce patient stress.

## High-quality Fundus Autofluorescence (FAF) Photography

Canon's autofluorescence (FAF) photography uses a high-sensitivity digital camera unit capable of taking realistic, high-quality images.

## Opacity Suppression

An image processing method that makes optic nerve heads and blood vessels easier to see in unclear fundus images with cloudy optic media, such as with cataracts.



## Panorama Image Creation \* Optional

The software automatically stitches together multiple images taken from any desired direction.



## CX-1 Specifications

Type	Mydriatic and Non-mydriatic
Photography mode	COLOR, RED FREE, COBALT, FLUO (fluorescein angiography), FAF (autofluorescence)
Retinal observation	Mydriatic mode: Optical viewfinder
	Non-mydriatic mode: LCD monitor mounted on rear of digital camera unit
Field angle	Mydriatic mode: 50°
	Non-mydriatic mode: 45°
Magnification	2x (digital)
Mounted digital camera	CX-1 digital camera unit
Patient diopter compensation range	Without compensation lens: -10 D to +15 D
	When using negative compensation lens: -31 D to -7 D
	When using positive compensation lens: +11 D to +33 D
Working distance	35 mm
Focusing	Split-line alignment type
Working distance adjustment	Working dots
Eye fixation lamp	Mydriatic mode: External eye fixation lamp
	Non-mydriatic mode: External eye fixation lamp, internal eye fixation lamp
Observation light source	Mydriatic mode: Halogen lamp
	Non-mydriatic mode: Infrared LED
Photography light source	Xenon tube
Operating range	Stage: Front/back: 65 mm, Left/right: 110 mm
	Main unit up/down: 90 mm
Panning range	Left/right: 30°
Tilting range	Up: 15°
	Down: 10°
Operating environment	Temperature: 10°C to 35°C
	Humidity: 30% to 80% RH
Dimensions	W320 mm x L531 mm x H577 mm
Weight	26 kg

## Main Components

CX-1 main unit
Digital camera
External eye fixation lamp
Digital camera cover
Objective lens cap
Chin rest paper (100 sheets)
Dust cover

## Options

Internal eye fixation target

\* Specifications and appearance are subject to change without notice.



### ■ CANON SINGAPORE PTE. LTD.

Medical Equipment Product Division

1 Fusionopolis Place, #14-10 Galaxis, Singapore 138522

Telephone: +65-6796-3549

<https://sg.canon/en/business/products>



General Safety Information

- Make sure you read the manual before using the instrument.
- Use the correct power source and voltage shown.
- Make sure the instrument is properly connected to ground. Failure to do so may result in an electric shock if there is a malfunction or ground leakage.

